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Distributore in Italia: **SIEM - NOVA SRL** - 20089 ROZZANO (MI) - VIA VARALLI 1 (ITA) - Italy

Produttore: Elmaslar İmalat Tibbi Cihazlar İthalat İhracat A.Ş. AĞAÇ METAL SAN. KOOP.1436 SK NO: 14 OSTİM MAH.ANKARA TURKEY

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Emission test	Compliance	Electromagnetic environment – guidance
Emissions RF CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Emissions RF CISPR 11	Class B	The device is suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions used for domestic purposes. IEC 61000-3-2	Class A	
Voltage fluctuations/ flickeremissions IEC 61000-3-3	Complies	

Immunity test	IEC 60601 Testlevel	Compliancelevel	Electromagnetic environment – guidance
Electrostatic discharges IEC 61000-4-2	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30
Electrical fast transient/burst IEC 61000-4-4	Power supplylines: ±2 kV Input/output lines: ±1 kV 100 kHz repetition frequency	Power supplylines: ±2 kV Input/output lines: ±1 kV 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surges Line-to-line IEC 61000-4-5	± 0,5 kV, ± 1 kV	± 0,5 kV, ± 1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surges Line-to-ground IEC 61000-4-5	± 0,5 kV, ± 1 kV, ± 2 kV	± 0,5 kV, ± 1 kV, ± 2 kV	
Voltage dips IEC 61000-4-11	0% UT; 0,5 cycle 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	0% UT; 0,5 cycle 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	Mains power quality should be that of a typical commercial or hospital environment.
	0% UT; 1 cycle and 70% UT; 25/30 cycle Single phase: 0°	0% UT; 1 cycle and 70% UT; 25/30 cycle Single phase: 0°	
Voltage interruptions IEC 61000-4-11	0% UT; 250/300 cycle	0% UT; 250/300 cycle	
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Not: UT is the a.c. mains voltage prior to application of the test level.

Conducted RF IEC 61000-4-6	150kHz to 80MHz: 3Vrms 6 Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	150kHz to 80MHz: 3Vrms 6 Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the device including cables, than the recommended separation distance calculated from the equation applicable to the frequency the transmitter. Recommended separation distance: $d = 0.35 \sqrt{P}$ $d = 1.2 \sqrt{80 \text{ MHz}}$ ile 800 MHz $d = 1.2 \sqrt{800 \text{ MHz}}$ ile 2.7 GHz $d = 2.3 \sqrt{P}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	10V/m 80% Am at 1kHz	10V/m 80% Am at 1kHz	

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Distributore in Italia: **SIEM - NOVA SRL** - 20089 ROZZANO (MI) - VIA VARALLI 1 (ITA) - Italy

Produttore: Elmaslar İmalat Tibbi Cihazlar İthalat İhracat A.Ş. AĞAÇ METAL SAN. KOOP.1436 SK NO: 14 OSTİM MAH.ANKARA TURKEY

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Note 1: 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

c. Calibration for current injection clamps shall be performed in a 150 Ω system.

d. If the frequency stepping skips over an ISM or amateur band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.

e. Applicable to ME EQUIPMENT and ME SYSTEMS with RATED input current less than or equal to 16 A / phase and ME EQUIPMENT and ME SYSTEMS with RATED input current greater than 16 A / phase.

Recommended separation distances between portable and mobile RF communications equipment and the device

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter(m)			
	150 kHz to 80 MHz (out ISM and amateur radio bands) $d= 1.2 \sqrt{P}$	150 kHz to 80 MHz (out ISM and amateur radio bands) $d= 0.6 \sqrt{P}$	80 MHz to 800 MHz $d= 1.2 \sqrt{P}$	800 MHz to 2.7 GHz $d= 2.3 \sqrt{P}$
0,01	0,12	0,06	0,12	0,23
0,1	0,38	0,19	0,38	0,73
1	1,2	0,6	1,2	2,3
10	3,8	1,9	3,8	7,3
100	12	6	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

Note 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Radiated RF	Test Frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	Immunity Test Level(V/m)
IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipments)	385	380-390	TETRA 400	Pulse modulation b) 18 Hz	1,8	0,3	27
	450	380-390	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0,3	28
	710 745 780	704-787	LTE BAND1 13, 17	Pulse modulation b) 217 Hz	0,2	0,3	9
	810 870 930	800-960	GSM 800/900, TETRA 800, IDEN 820, LTE Bant 5	Pulse modulation b) 18 Hz	2	0,3	28
	1720 1845 1970	1700-1990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Bant1,3,4,25, UMTS	Pulse modulation b) 217 Hz	2	0,3	28
	2450	2400-2570	Bluetooth, WLAN 802.11b/g/n, RFID 2450, LTE Bant 7	Pulse modulation b) 217 Hz	2	0,3	28
	5240 5240 5785	5100-5800	WLAN 802.11a/n	Pulse modulation b) 217 Hz	0,2	0,3	9

Note: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

a. For some services, only the uplink frequencies are included.

b. The carrier shall be modulated using a 50 % duty cycle square wave signal.

c. As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation:

$$E= 6 / d \sqrt{P}$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.